


# MAXICORE-M® BIOPSY GUN NEEDLE

## MAXICORE-M® BİYOPSİ TABANCA İĞNESİ

- ↳ Echogenic needle tip
- ↳ 17 mm sample notch
- ↳ Depth marked canula
- ↳ Color code
- ↳ Compatible with Maxicore-M® Reusable Biopsy Gun and Bard Magnum® Reusable Gun (\*)

(\*) Bard Magnum® is a Trade Mark of Bard Inc.



- ↳ Ekojenik uç
- ↳ 17mm örnek alma haznesi
- ↳ Derinlik işaretli kanül
- ↳ Renk kodu
- ↳ Maxicore-M™ ve Bard Magnum™ Reusable Biyopsi Tabancası ile uyumludur

Bard Magnum® is a Trade Mark of Bard Inc.

Ref. No	Diameter (G)	Length (cm)
GPN1410-M	14	10
GPN1415-M	14	15
GPN1420-M	14	20
GPN1610-M	16	10
GPN1615-M	16	15
GPN1620-M	16	20
GPN1810-M	18	10
GPN1815-M	18	15
GPN1820-M	18	20
GPN1825-M	18	25
GPN1830-M	18	30